

Product data sheet

Medical Face Mask


REF:	204
Product:	Medical Face Mask
Manufacturer:	DACH Schutzbekleidung GmbH & Co. KG



Categorization:	Class I Medical Device Regulation (EU) 2017/745; UMDNS Code: 12-458. Type IIR according to European standard EN 14683:2019 for surgical masks		
Product description:	Ergonomically formed face mask with easily adjustable nose piece. The outer and the inner layer consist of physiologically harmless soft non-woven fabric. Joined by ultrasonic welding and equipped with elastic ear loops.		
Material:	Outer layer: Polypropylene Inner layer: Soft Polypropylene, skin friendly, allergy free Filter media: Polypropylene Ear loop: soft polyester and spandex Nose piece: Plastic coated iron wire		
Dermatological tolerance:	No intolerance reaction.		
Application:	All used materials are irritant free and odorless, qualified for medical use. By wearing this mask, the user protects his environment.		
Color:	Black		
Packaging unit:	Units	Packaging	GTIN
	50 pcs.	Box	4049825006459
	2000 pcs.	Shipping carton	
Storage:	Store in a dry environment out of direct sunlight in original packaging, see also packaging. Shelf life: 5 years		



DACH Schutzbekleidung
GmbH & Co. KG

Disposal:	<p>This product can be thermally recycled or disposed in a controlled landfill without releasing toxic substances.</p> <p>The disposal of contaminated products must be determined by its contamination. Please respect according laws and regulations.</p>
Indication:	<p>Non sterile. </p>